

Tecan M Nano: Fluorescence + OD600

platereader
scan for the
full protocol



PROTOCOL

- 1. Fill the plate.** P200 multichannel, **100 μ L** per well into a **black-walled, clear-bottom 96-well plate**. Each culture **2 $\times$** . Skip any culture that did not grow.
 - Tips all seated, pipette up and down to resuspend, no air gaps.
Volume sets the signal — pipette accurately.
- 2. Start up, in this order, from everything off:** the **computer** (PIN from a supervisor), then the **instrument** (toggle **on the back**), then **iControl**.
- 3. Connect:** choose **Infinite 200 Pro**, OK. It always takes a while.
 - **Will not connect?** End iControl (Task Manager if it hangs), switch everything off, and start again **in the same order**.
- 4. Default script.** Select only the **wells you filled**.
- 5. Add Absorbance: 600 nm.** Add **Fluorescence Intensity:** ex **483 nm**, em **525 nm**, gain **40**. **Leave everything else at its default.**
- 6. Open the tray** (iControl or the corner button). Plate **in the slot, A1 top left**. Close the same way. **Never push or pull the door by hand.**
- 7. Start.** Results open in **Excel**, absorbance and fluorescence on separate sheets. **File** $\rightarrow$ **Save As** to the desktop with a real name, then USB stick or email.
- 8. Plate out first**, close the tray. Quit Excel; quit iControl — **yes** disconnect, **no** save script. Shut down the computer, **then** the instrument.
- 9. Bleach and rinse the plate.**

Reads OVER? Out of range — lower the gain and read again. Numbers weak — raise it.