
PROTOCOL

1. For each clone, **top-label a 1.5 mL tube** with the clone ID from the labsheet. This is the tube you will elute into.
 2. **Pellet 4 mL** of saturated culture. Spin 1 min, pour off the supernatant — **bleach it** before it goes down the drain.
 3. **Resuspend** in **250 µL P1**. No clumps. The RNase A must already be in the P1.
 4. **Lyse** with **250 µL P2**. Mix gently — **do not vortex**. It goes clear and viscous.
 5. **Neutralise** with **350 µL N3**. Invert thoroughly. A white precipitate forms.
 6. **Spin 5 min** at max speed to pellet the debris.
 7. Transfer the **supernatant** to a blue QIAprep column. Spin **15 s**.
 8. Add **500 µL PB**, spin. Add **750 µL PE**, spin.
 9. Discard the flow-through and **spin 90 s to dry**. PE is 70% ethanol and carryover ruins everything downstream.
 10. Move the column to your labelled tube. Add **50 µL EB** to the **centre of the membrane** and spin **45 s**.
- Label **the top and the side** of each tube — a tube marked only on the lid goes anonymous the moment someone puts the lids down.
 - **A clone that yields no DNA was probably a satellite colony**. Note it and exclude it.