

---

## PROTOCOL

---

1. For each cleanup, **top-label a 1.5 mL tube** with the name indicated on the labsheet. This is the tube you will elute into.
  2. Add **180 µL ADB** (brown bottle) to the reaction and mix.
    - Fragments under **250 bp**: use **1 part ADB + 3 parts isopropanol** instead of ADB alone, or the fragment washes straight through.
  3. Transfer to a **Zymo column** in a collection tube. Spin **15 s** at full speed, discard the flow-through.
  4. Add **200 µL PE**. Spin **15 s**, discard the flow-through.
  5. Add **200 µL PE** again. Spin **15 s**, discard the flow-through.
  6. Spin **90 s** at full speed to **dry the column**. PE is 70% ethanol and any carryover inhibits the ligase.
  7. Move the column to your labelled tube. Add **25 µL EB** slowly to the **centre of the membrane** — do not let it run down the walls.
  8. Spin **45 s** to elute. Discard the column.
- **Know where your DNA is at every step.** Before each spin, say which half you are keeping.