

## REACTION — 10 $\mu$ L

**Label the tube first**, before any liquid goes in. The **top label** is the number from your labsheet.

**6** ddH<sub>2</sub>O (white rack)

**1** 10 $\times$  T4 DNA ligase buffer (red)

**1** each DNA fragment (2 of them, 2  $\mu$ L total)

**0.5** T4 DNA ligase

**0.5** Type IIS enzyme — from your labsheet

**10  $\mu$ L total**

Top to bottom. **Enzymes last.** Mix well, quick spin.

## ENZYME

**Not always BsaI.** It must match the sites designed into your fragments. Take it from your labsheet, not the nearest tube.

## RUN

All reactions from your section in **one block**. Run **GG1**.

Carry the whole 10  $\mu$ L into the transformation. **Do not split it.**