

WHAT THIS IS FOR

Checking that a PCR produced a product and that it is the right size. You are **not purifying** — nothing is recovered from this gel.

Slabs are **1% agarose in 1× TAE**, cut from the pre-made gels in the fridge.

PREPARE EACH SAMPLE

8 loading dye (tube marked *load*)

3 PCR product

Mix, quick spin.

Marker, one per section: 8 µL loading dye into a tube of marker (yellow), mix, spin.

SET UP AND LOAD

1. Cut a slab with enough wells for every sample, the marker, and a few spare.
2. Put it in the rig. Fill with **1× TAE** until it is just covered.
3. Brace the gel with the plastic wedge so it cannot float.
4. Write the sample order on a paper strip; line the tubes up to match.
5. Load **9 µL** per well with a P20.

RUN

Run to red. DNA is negative and moves to the red electrode, so load at the black end.

Run at **175 V** for **~10 min**, or until the blue front is **2/3–3/4 down the gel**.

Watch the gel, not the clock. 175 V for 10 min has not been re-measured on the B144 power supply and rigs. The dye front is the real endpoint. Tell your supervisor what actually worked.

IMAGE

1. Move the gel to the imager and lay the label strip above it.
2. Photograph it **with your phone through the orange filter**.
3. Name the file by date and time, then upload it.

READING IT

- Smaller fragments move faster.
- A **library** gives a smear, not a sharp band — its members differ in length.
- Loading dye does two jobs: glycerol sinks the sample into the well, and the colour tracks the run. It does not bind DNA.
- Load below the full volume on purpose — 9 µL of 11 µL — so you do not draw an air gap.
- Do not puncture a well. The sample runs out the bottom; use another lane.